
# You have the hands of a surgeon

a surgeon and the eyes of an artist.

That's why your OVD is HEALON®

Trust the HEALON® OVDs to support your craft and provide premium protection for patient eyes throughout the cataract surgery









#### **SPEC SHEET**

## Healon GV® PRO OVD

### Ophthalmic Viscosurgical Devices (OVD)

#### **Specifications**



	Healon GV® PRO OVD
Classification	Cohesive
Composition	1.8% Sodium hyaluronate¹
SKUs	10270014 - (CE Central EN, ES, PT, EL) 10290014 - (CE Central IT, DE, NL, FR) 10291014 - (CE Nordics DK, FI, SE, NO) 10272014 - (CE East CS, HR, SR, SK) 10282014 - (CE East HU, RU) 10292014 - (CE East BG, SL, RO, TR) 10278014 - (CE Arabic)
Volume	0.85 ml¹
Zero Shear Viscosity	2.60 - 3.50 Log (Pas) <sup>2</sup>
Average Molecular Weight	3,200,000 Daltons <sup>1</sup>
Storage Temperature	+2°C to +8°C1

All **HEALON®** OVDs are latex free.

#### References

- 1. **HEALON GV PRO** OVD DfU Doc. 52-0480-00. 2020-04-29. REF2020CT4146.
- 2. Evaluation of the rheological properties of *HEALON GV PRO OVD SEY 0389-02*. REF2019CT4044.

Only for healthcare professionals. Please read the Directions for Use for Important Safety Information and consult our specialists if you have any questions.