
You have the hands of a surgeon

a surgeon and the eyes of an artist.

That's why your OVD is HEALON®

Trust the HEALON® OVDs to support your craft and provide premium protection for patient eyes throughout the cataract surgery





FROM THE MAKERS OF

TECNIS° Family of IOLs

Ophthalmic Viscosurgical Devices (OVD)

Specifications



2.3% Sodium hyaluronate

| | Healon5® PRO OVD |
|--------------------------|--|
| Classification | Viscoadaptive |
| Composition | 2.3% Sodium hyaluronate ¹ |
| SKUs | 10290015 (CE Central) 10291015 (CE Nordics) 10292015 (CE East) |
| Volume | 0.60 ml¹ |
| Zero Shear Viscosity | 1,300,000 - 12,600,000 mPas² |
| Average Molecular Weight | 3,200,000 Daltons ¹ |
| Storage Temperature | +2°C to +8°C1 |

All **HEALON®** OVDs are latex free.

References

1. HEALON5 PRO OVD DfU - CE Central - Doc. 52-469-00 - Rev 2020-01-14. REF2020CT4204

Only for healthcare professionals. Please read the Directions for Use for Important Safety Information and consult our specialists if you have any questions.

^{2.} D0F2016CT0020 - Evaluation of rheological properties of *HEALON* PRO OVD and *HEALON5* PRO OVD. 2016-06-13.